




Product Specification
Transfer pipette, 6 ml, (LxW): 152 x 15 mm, LD-PE,
transparent

Page 1



Product description

Order number	86.1174.300
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 152 x 15 mm, material: LD-PE, transparent, droplet size: 40 µl - 55 µl

Product characteristics

Graduation	no
------------	----

Size

Sample volume	6 ml
Width of product	15 mm
Length of product	152 mm
Drop size	40 µl - 55 µl

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Transfer pipette, 6 ml, (LxW): 152 x 15 mm, LD-PE,
transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	400
Piece(s) / inner box	400
Piece (s) / outer case	4000
Piece(s) / pallet	44800
Case volume	0.0994 cbm
Weight of product	0.002 kg
Weight of case	9.07 kg
EAN of inner box	4038917378505
EAN case	4038917378499

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

