




**Product Specification**  
**Transfer pipette, 6 ml, (LxW): 146 x 15 mm, LD-PE,**  
**transparent**

**Page 1**



**Product description**

Order number	86.1175.010
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 146 x 15 mm, material: LD-PE, transparent, droplet size: 20 µl - 30 µl, sterile, 10 piece(s)/blister

**Product characteristics**

Graduation	no
------------	----

**Size**

Sample volume	6 ml
Width of product	15 mm
Length of product	146 mm
Drop size	20 µl - 30 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	blister
Piece(s) / inner box	250
Piece (s) / outer case	1000
Piece(s) / pallet	28000
Depth of box	261 mm
Width of box	161 mm
Height of box	231 mm
Depth of case	533 mm
Width of case	333 mm
Height of case	250 mm
Case volume	0.0444 cbm
Weight of product	0.0015 kg
Weight of case	2.71 kg
EAN of inner box	4038917353168
EAN case	4038917353137

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

