



**Product Specification**  
**Transfer pipette, 1 ml, (LxW): 87 x 10 mm, LD-PE,**  
**transparent**

**Page 1**



**Product description**

Order number	86.1179.300
Product description	Transfer pipette, nominal volume: 1 ml, (LxW): 87 x 10 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl, 500 piece(s)/case

**Product characteristics**

Graduation	no
------------	----

**Size**

Sample volume	1 ml
Width of product	10 mm
Length of product	87 mm
Drop size	30 µl - 45 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 1 ml, (LxW): 87 x 10 mm, LD-PE,**  
**transparent**

**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	500
Piece (s) / outer case	4000
Piece(s) / pallet	234000
Case volume	0.041 cbm
Weight of product	0.001 kg
Weight of case	4.55 kg
EAN of inner box	4038917378598
EAN case	4038917378581

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

