



**Product Specification**  
**Sarpette® M, single-channel, variable, 100 µl - 1,000 µl, blue, 1 piece(s)/case**

**Page 1**



**Product description**

Order number	90.3100.000
Product description	Sarpette® M, single-channel, variable, 100 µl - 1,000 µl, blue, 1 piece(s)/case

**Product characteristics**

Volume of work	100 µl - 1.000 µl
Coding colour	blue
Autoclavable (max. °C)	121 °C

**Material & colours**

Colour of product	grey
-------------------	------

**Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Sarpette® M, single-channel, variable, 100 µl -**  
**1,000 µl, blue, 1 piece(s)/case**

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Depth of box	315 mm
Width of box	90 mm
Height of box	30 mm
Weight of product	0.088 kg
EAN of inner box	4038917937061
EAN case	4038917937078

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

