



### **Product description**

Order number	90.3112.010
Product description	Sarpette® M, 12-channel, variable, 0.5 µl - 10 µl, grey, 1 piece(s)/case

### **Product characteristics**

Volume of work	0,5 µl - 10 µl
Coding colour	grey
Autoclavable (max. °C)	121 °C

### **Material & colours**

Colour of product	grey
-------------------	------

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Depth of box	288 mm
Width of box	150 mm
Height of box	28 mm
Weight of product	0.185 kg
EAN of inner box	4038917937610
EAN case	4038917937627

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

