



## Product description

Order number	70.1181
Product description	Pipette tip, working volume: 1,000 µl, transparent, suitable for MLA pipettes, 250 piece(s)/bag

## Product characteristics

Graduation	no
Suitable for	MLA pipettes
Volume of work	1.000 µl
Filter	no
Autoclavable (max. °C)	121 °C

## Size

Max. volume	1000 µl
Length of product	86 mm

## Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

## Purity & certification

Product category	no medical device   IVD
Batched	no

## Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	5000
Piece(s) / pallet	180000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0007 kg
Weight of case	4.41 kg
EAN of inner box	4038917085656

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Pipette tip, 1,000 µl, transparent, 250 piece(s)/bag**

**Page 2**

EAN case

4038917004299

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**