




**Product Specification**  
**Filter tip, 1,250 µl, transparent, Biosphere® , 96**  
**piece(s)/box**

**Page 1**



**Product description**

Order number	70.1186.210
Product description	Filter tip, working volume: 1,250 µl, transparent, filling level rings, Biosphere® , suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 96 piece(s)/box

**Product characteristics**

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design
Volume of work	1.250 µl
Filter	yes
Coding colour	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Filter tip, 1,250 µl, transparent, Biosphere® , 96  
piece(s)/box**

**Page 2**

**Size**

Max. volume	1250 µl
Length of product	103 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Filter material	Polyethylene (PE)

**Purity & certification**

Product category	no medical device   IVD
Purity standard	Biosphere®
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	768
Type of smallest subpackaging	box
Piece(s) / inner box	96
Piece (s) / outer case	768
Piece(s) / pallet	90000
Depth of case	270 mm
Width of case	121 mm
Height of case	413 mm
Case volume	0.0135 cbm
Weight of product	0.0022 kg
Weight of case	2.08 kg
EAN of inner box	4038917126496
EAN case	4038917080354

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-16

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**