



Product description

Order number	70.1186.410
Product description	Filter tip, working volume: 1,250 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 384 piece(s)/bag

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design
Volume of work	1.250 µl
Filter	yes
Coding colour	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Filter tip, 1,250 µl, transparent, PCR Performance
Tested, 384 piece(s)/bag

Page 2

Size

Max. volume	1250 µl
Length of product	103 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Filter material	Polyethylene (PE)

Purity & certification

Product category	no medical device IVD
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	3840
Type of smallest subpackaging	bag
Piece(s) / inner box	384
Piece (s) / outer case	3840
Piece(s) / pallet	134400
Depth of case	343 mm
Width of case	248 mm
Height of case	438 mm
Case volume	0.0373 cbm
Weight of product	0.0011 kg
Weight of case	4.82 kg
EAN of inner box	4038917126502
EAN case	4038917080361

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT