



### **Product description**

Order number	70.1187.002
Product description	Pipette tip, working volume: 10 ml, blue, filling level rings, suitable for Sarpette® M, 100 piece(s)/bag

### **Product characteristics**

Graduation	yes
Suitable for	Sarpette® M
Volume of work	10 ml
Filter	no
Autoclavable (max. °C)	121 °C

### **Size**

Max. volume	10 ml
Length of product	144 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Pipette tip, 10 ml, blue, 100 piece(s)/bag**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	blue

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	20000
Depth of case	457 mm
Width of case	381 mm
Height of case	269 mm
Case volume	0.0468 cbm
Weight of product	0.0048 kg
Weight of case	5.4 kg
EAN of inner box	4038917951531
EAN case	4038917951548

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-05

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**