




Product Specification
Pipette tip, extra long, 200 µl, transparent, 96
piece(s)/box

Page 1



Product description

Order number	70.1189.105
Product description	Pipette tip, working volume: 200 µl, transparent, filling level rings, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and products of identical design, extra long, 96 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and products of identical design
Volume of work	200 µl
Filter	no
Coding colour	yellow
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Pipette tip, extra long, 200 µl, transparent, 96
piece(s)/box**

Page 2

Size

Max. volume	200 µl
Length of product	77 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	960
Type of smallest subpackaging	box
Piece(s) / inner box	96
Piece (s) / outer case	960
Piece(s) / pallet	120960
Depth of case	263 mm
Width of case	121 mm
Height of case	413 mm
Case volume	0.0131 cbm
Weight of product	0.0015 kg
Weight of case	1.94 kg
EAN of inner box	4038917314497
EAN case	4038917314480

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-03

This document was prepared by EDP
support and is valid without signature.



SARSTEDT