




Product Specification
Filter tip, 200 µl, transparent, Biosphere® , 96
piece(s)/box

Page 1



Product description

Order number	70.1189.215
Product description	Filter tip, working volume: 200 µl, transparent, filling level rings, Biosphere® , suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Finnpiquette and products of identical design, 96 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson and Finnpiquette and products of identical design
Volume of work	200 µl
Filter	yes
Coding colour	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Filter tip, 200 µl, transparent, Biosphere® , 96
piece(s)/box**

Page 2

Size

Max. volume	200 µl
Length of product	77 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Filter material	Polyethylene (PE)

Purity & certification

Product category	no medical device IVD
Purity standard	Biosphere®
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	960
Type of smallest subpackaging	box
Piece(s) / inner box	96
Piece (s) / outer case	960
Piece(s) / pallet	120960
Depth of case	241 mm
Width of case	89 mm
Height of case	406 mm
Case volume	0.0087 cbm
Weight of product	0.0014 kg
Weight of case	1.94 kg
EAN of inner box	4038917314589
EAN case	4038917314572

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

