



Product description

Order number	70.1193.151
Product description	Filter tip, working volume: 1,000 µl, black, suitable for application on SARSTEDT IVARO AP devices, 480 piece(s)/blister

Product characteristics

Suitable for	application on SARSTEDT IVARO AP devices
Volume of work	1.000 µl
Filter	yes

Size

Max. volume	1000 µl
-------------	---------

Material & colours

Product material	Polypropylene (PP)
Colour of product	black

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	3840
Type of smallest subpackaging	blister
Piece(s) / inner box	480
Piece (s) / outer case	3840
Piece(s) / pallet	61440
Depth of case	500 mm
Width of case	335 mm
Height of case	200 mm
Case volume	0.0754 cbm
Weight of product	0.002 kg
Weight of case	7.95 kg
EAN of inner box	4038917913447
EAN case	4038917913454

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

