



Product description

| | |
|---------------------|--|
| Order number | 70.3010.100 |
| Product description | Pipette tip, working volume: 10 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 480 piece(s)/StackPack |

Product characteristics

| | |
|----------------|---|
| Graduation | yes |
| Suitable for | SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design |
| Volume of work | 10 µl |
| Filter | no |
| Coding colour | grey |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|-------|
| Max. volume | 10 µl |
| Length of product | 31 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1920 |
| Type of smallest subpackaging | StackPack |
| Piece(s) / inner box | 1920 |
| Piece (s) / outer case | 1920 |
| Piece(s) / pallet | 307200 |
| Depth of case | 290 mm |
| Width of case | 143 mm |
| Height of case | 225 mm |
| Case volume | 0.0093 cbm |
| Weight of product | 0.0004 kg |
| Weight of case | 1.69 kg |
| EAN of inner box | 4038917894432 |
| EAN case | 4038917894425 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

