



### Product description

Order number	70.3010.305
Product description	Pipette tip, working volume: 10 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 piece(s)/SingleRefill

### Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design
Volume of work	10 µl
Filter	no
Coding colour	grey
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Pipette tip, 10 µl, transparent, Biosphere® plus,**  
**96 piece(s)/SingleRefill**

**Page 2**

**Size**

Max. volume	10 µl
Length of product	31 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Batched	yes

**Packaging**

Minimum order qty.	1920
Type of smallest subpackaging	SingleRefill
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	115200
Depth of box	291 mm
Width of box	105 mm
Height of box	148 mm
Depth of case	437 mm
Width of case	307 mm
Height of case	170 mm
Case volume	0.0228 cbm
Weight of product	0.0006 kg
Weight of case	5.16 kg
EAN of inner box	4038917894951
EAN case	4038917894968

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-16

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**