



Product description

Order number	70.3010.365
Product description	Pipette tip, working volume: 2.5 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 piece(s)/SingleRefill

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design
Volume of work	2,5 µl
Filter	yes
Coding colour	grey

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Pipette tip, 2.5 µl, transparent, Biosphere® plus,
96 piece(s)/SingleRefill

Page 2

Size

Max. volume	2.5 µl
Length of product	31 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Filter material	Polyethylene (PE)

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order qty.	1920
Type of smallest subpackaging	SingleRefill
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	115200
Depth of box	530 mm
Width of box	111 mm
Height of box	148 mm
Depth of case	547 mm
Width of case	325 mm
Height of case	245 mm
Case volume	0.0436 cbm
Weight of product	0.0006 kg
Weight of case	5.16 kg
EAN of inner box	4038917895040
EAN case	4038917895057

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT