



**Product description**

Order number	70.3020
Product description	Pipette tip, working volume: 20 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 1,000 piece(s)/bag

**Product characteristics**

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design
Volume of work	20 µl
Filter	no
Coding colour	grey
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Pipette tip, 20 µl, transparent, PCR Performance**  
**Tested, 1,000 piece(s)/bag**

**Page 2**

**Size**

Max. volume	20 µl
Length of product	46 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

**Packaging**

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	720000
Depth of case	533 mm
Width of case	263 mm
Height of case	140 mm
Case volume	0.0196 cbm
Weight of product	0.0002 kg
Weight of case	2.19 kg
EAN of inner box	4038917525862
EAN case	4038917525855

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-24

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**