



### Product description

Order number	70.3030.275
Product description	Filter tip, working volume: 100 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/box

### Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design
Volume of work	100 µl
Filter	yes
Coding colour	yellow
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Filter tip, 100 µl, transparent, Biosphere® plus,**  
**Low retention, 96 piece(s)/box**

**Page 2**

**Size**

Max. volume	100 µl
Length of product	51 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Filter material	Polyethylene (PE)

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Batched	yes

**Packaging**

Minimum order qty.	1920
Type of smallest subpackaging	box
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	96000
Depth of box	330 mm
Width of box	111 mm
Height of box	148 mm
Depth of case	348 mm
Width of case	325 mm
Height of case	245 mm
Case volume	0.0277 cbm
Weight of product	0.0011 kg
Weight of case	5.16 kg
EAN of inner box	4038917846165
EAN case	4038917846172

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-25

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**