



**Product description**

Order number	70.3030.365
Product description	Filter tip, working volume: 20 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/SingleRefill

**Product characteristics**

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design
Volume of work	20 µl
Filter	yes
Coding colour	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Max. volume	20 µl
Length of product	51 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Filter material	Polyethylene (PE)

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Batched	yes

### **Packaging**

Minimum order qty.	1920
Type of smallest subpackaging	SingleRefill
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	115200
Depth of box	530 mm
Width of box	111 mm
Height of box	148 mm
Depth of case	547 mm
Width of case	325 mm
Height of case	245 mm
Case volume	0.0436 cbm
Weight of product	0.0008 kg
Weight of case	5.16 kg
EAN of inner box	4038917895002
EAN case	4038917895019

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

