



Product description

Order number	70.3031.110
Product description	Pipette tip, working volume: 250 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 480 piece(s)/StackPack

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design
Volume of work	250 µl
Filter	no
Coding colour	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Max. volume	250 µl
Length of product	51 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	1920
Type of smallest subpackaging	StackPack
Piece(s) / inner box	1920
Piece (s) / outer case	1920
Piece(s) / pallet	291840
Depth of case	290 mm
Width of case	157 mm
Height of case	225 mm
Case volume	0.0102 cbm
Weight of product	0.0003 kg
Weight of case	1.69 kg
EAN of inner box	4038917845533
EAN case	4038917845540

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

