



Product description

| | |
|---------------------|---|
| Order number | 70.3031 |
| Product description | Pipette tip, working volume: 250 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 1,000 piece(s)/bag |

Product characteristics

| | |
|------------------------|---|
| Graduation | yes |
| Suitable for | SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design |
| Volume of work | 250 µl |
| Filter | no |
| Coding colour | yellow |
| Autoclavable (max. °C) | 121 °C |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|--------|
| Max. volume | 250 µl |
| Length of product | 51 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 10000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 10000 |
| Piece(s) / pallet | 480000 |
| Depth of case | 388 mm |
| Width of case | 378 mm |
| Height of case | 215 mm |
| Case volume | 0.0315 cbm |
| Weight of product | 0.0003 kg |
| Weight of case | 3.39 kg |
| EAN of inner box | 4038917526029 |
| EAN case | 4038917526012 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

