



Product description

| | |
|---------------------|---|
| Order number | 70.3040.100 |
| Product description | Pipette tip, working volume: 300 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 480 piece(s)/StackPack |

Product characteristics

| | |
|----------------|---|
| Graduation | yes |
| Suitable for | SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design |
| Volume of work | 300 µl |
| Filter | no |
| Coding colour | orange |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Pipette tip, 300 µl, transparent, PCR Performance
Tested, 480 piece(s)/StackPack

Page 2

Size

| | |
|-------------------|--------|
| Max. volume | 300 µl |
| Diameter | 7 mm |
| Length of product | 60 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |

Purity & certification

| | |
|---------------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| Certificate of conformity | yes |
| CE certified | CE |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1920 |
| Type of smallest subpackaging | StackPack |
| Piece(s) / inner box | 1920 |
| Piece (s) / outer case | 1920 |
| Piece(s) / pallet | 245760 |
| Depth of case | 290 mm |
| Width of case | 186 mm |
| Height of case | 225 mm |
| Case volume | 0.0121 cbm |
| Weight of product | 0.0007 kg |
| Weight of case | 1.69 kg |
| EAN of inner box | 4038917433433 |
| EAN case | 4038917433426 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT