



Product description

Order number	70.3040.375
Product description	Filter tip, working volume: 300 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 96 piece(s)/SingleRefill

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design
Volume of work	300 µl
Filter	yes
Coding colour	orange

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Filter tip, 300 µl, transparent, Biosphere® plus,
Low retention, 96 piece(s)/SingleRefill**

Page 2

Size

Max. volume	300 µl
Diameter	7 mm
Length of product	60 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
Certificate of conformity	yes
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	1920
Type of smallest subpackaging	SingleRefill
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	103680
Depth of box	333 mm
Width of box	105 mm
Height of box	148 mm
Depth of case	437 mm
Width of case	351 mm
Height of case	172 mm
Case volume	0.0264 cbm
Weight of product	0.0009 kg
Weight of case	2.81 kg
EAN of inner box	4038917433167
EAN case	4038917433150

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT