




## Product Specification

**Pipette tip, 1,000 µl, transparent, PCR  
Performance Tested, Low retention, 480  
piece(s)/StackPack**

**Page 1**



### Product description

Order number	70.3050.110
Product description	Pipette tip, working volume: 1,000 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiquette, Biohit and Brand and products of identical design, 480 piece(s)/StackPack

### Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiquette, Biohit and Brand and products of identical design
Volume of work	1.000 µl
Filter	no
Coding colour	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Pipette tip, 1,000 µl, transparent, PCR  
Performance Tested, Low retention, 480  
piece(s)/StackPack**

**Page 2**

### Size

Max. volume	1000 µl
Diameter	9 mm
Length of product	72 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Product category	In vitrodiagnostic device, CE
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

### Packaging

Minimum order qty.	1920
Type of smallest subpackaging	StackPack
Piece(s) / inner box	1920
Piece (s) / outer case	1920
Piece(s) / pallet	149760
Depth of case	303 mm
Width of case	295 mm
Height of case	225 mm
Case volume	0.0201 cbm
Weight of product	0.0012 kg
Weight of case	2.79 kg
EAN of inner box	4038917424813
EAN case	4038917424806

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-26

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**