




**Product Specification**  
**Pipette tip, 1,000 µl, transparent, PCR**  
**Performance Tested, 96 piece(s)/box**

**Page 1**



**Product description**

Order number	70.3050.200
Product description	Pipette tip, working volume: 1,000 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/box

**Product characteristics**

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design
Volume of work	1.000 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Pipette tip, 1,000 µl, transparent, PCR**  
**Performance Tested, 96 piece(s)/box**

**Page 2**

**Size**

Max. volume	1000 µl
Diameter	9 mm
Length of product	72 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

**Packaging**

Minimum order qty.	1920
Type of smallest subpackaging	box
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	61440
Depth of box	530 mm
Width of box	111 mm
Height of box	148 mm
Depth of case	547 mm
Width of case	325 mm
Height of case	245 mm
Case volume	0.0436 cbm
Weight of product	0.0018 kg
Weight of case	4.98 kg
EAN of inner box	4038917424783
EAN case	4038917424776

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**