




Product Specification
Pipette tip, 1,000 µl, transparent, Biosphere®
plus, 96 piece(s)/box

Page 1



Product description

Order number	70.3050.205
Product description	Pipette tip, working volume: 1,000 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design
Volume of work	1.000 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Pipette tip, 1,000 µl, transparent, Biosphere®
plus, 96 piece(s)/box

Page 2

Size

Max. volume	1000 µl
Diameter	9 mm
Length of product	72 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In vitro diagnostic device, CE
Certificate of conformity	yes
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	1920
Type of smallest subpackaging	box
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	61440
Depth of box	530 mm
Width of box	111 mm
Height of box	148 mm
Depth of case	547 mm
Width of case	325 mm
Height of case	245 mm
Case volume	0.0436 cbm
Weight of product	0.0018 kg
Weight of case	4.99 kg
EAN of inner box	4038917424752
EAN case	4038917424745

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT