




Product Specification

Pipette tip, 1,250 µl, transparent, PCR
Performance Tested, Low retention, 480
piece(s)/StackPack

Page 1



Product description

Order number	70.3060.110
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiquette, Biohit and Brand and products of identical design, 480 piece(s)/StackPack

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiquette, Biohit and Brand and products of identical design
Volume of work	1.250 µl
Filter	no
Coding colour	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-24

This document was prepared by EDP
support and is valid without signature.



SARSTEDT



Product Specification

Pipette tip, 1,250 µl, transparent, PCR
Performance Tested, Low retention, 480
piece(s)/StackPack

Page 2

Size

Max. volume	1250 µl
Diameter	9 mm
Length of product	90 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	1920
Type of smallest subpackaging	StackPack
Piece(s) / inner box	1920
Piece (s) / outer case	1920
Piece(s) / pallet	184320
Depth of case	290 mm
Width of case	267 mm
Height of case	225 mm
Case volume	0.0174 cbm
Weight of product	0.0013 kg
Weight of case	2.97 kg
EAN of inner box	4038917432238
EAN case	4038917432221

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-24

This document was prepared by EDP support and is valid without signature.



SARSTEDT