




**Product Specification**  
**Filter tip, XL, 1,000 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/box**

**Page 1**



**Product description**

Order number	70.3060.275
Product description	Filter tip, working volume: 1,000 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, XL, 96 piece(s)/box

**Product characteristics**

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design
Volume of work	1.000 µl
Filter	yes
Coding colour	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Filter tip, XL, 1,000 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/box**

**Page 2**

**Size**

Max. volume	1000 µl
Diameter	9 mm
Length of product	90 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
Certificate of conformity	yes
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	1920
Type of smallest subpackaging	box
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	61440
Depth of box	530 mm
Width of box	111 mm
Height of box	148 mm
Depth of case	547 mm
Width of case	325 mm
Height of case	245 mm
Case volume	0.0436 cbm
Weight of product	0.002 kg
Weight of case	5.41 kg
EAN of inner box	4038917432290
EAN case	4038917432283

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-25

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**