



**Product Specification**  
**Replacement filter, for Automatic-Sarpette® , pore**  
**size: 0.2 µm**

**Page 1**



**Product description**

Order number	90.189.221
Product description	Replacement filter, for Automatic-Sarpette® , pore size: 0.2 µm, 5 piece(s)/bag

**Product characteristics**

Suitable for	Automatic-Sarpette®
Pore size	0,2

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	5
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	500
Piece(s) / pallet	54000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Replacement filter, for Automatic-Sarpette® , pore**  
**size: 0.2 µm**

Depth of case	308 mm
Width of case	223 mm
Height of case	195 mm
Case volume	0.0134 cbm
Weight of product	0.002 kg
Weight of case	1.21 kg
EAN of inner box	4038917419604
EAN case	4038917419598

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

