



### **Product description**

Order number	90.189.222
Product description	Silicon adapter, for Automatic-Sarpette® , translucent, silicone, 1 piece(s)/bag

### **Product characteristics**

Suitable for	Automatic-Sarpette®
--------------	---------------------

### **Material & colours**

Product material	silicone
Colour of product	translucent

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	100
Piece(s) / pallet	32400
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.014 kg
Weight of case	1.65 kg
EAN of inner box	4038917419086
EAN case	4038917419079

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

