




**Product Specification**  
**Serological pipette, with tip, plugged, 2 ml, sterile,**  
**25 piece(s)/bag**

**Page 1**



**Product description**

Order number	86.1252.025
Product description	Serological pipette, with tip, plugged, 2 ml, coding colour: green, material: PS, sterile, 25 piece(s)/bag

**Product characteristics**

Colour of print/label	green
Graduation	yes
Design	with tip, plugged
Volume of work	2 ml
Coding colour	green

**Size**

Diameter	6 mm
Length of product	270 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	1000
Piece(s) / pallet	81000
Depth of case	413 mm
Width of case	250 mm
Height of case	195 mm
Case volume	0.0201 cbm
Weight of product	0.0036 kg
Weight of case	4.21 kg
EAN of inner box	4038917127219
EAN case	4038917021098

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

