



Product description

Order number	65.716.300
Product description	Screw cap, natural, suitable for screw cap micro tubes, 500 piece(s)/bag

Product characteristics

Cap	flat cap
Closure type	screw cap
Suitable for	screw cap micro tubes

Material & colours

Product material	Polypropylene (PP)
Colour of product	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	10000
Piece(s) / pallet	720000
Depth of case	438 mm
Width of case	268 mm
Height of case	155 mm
Case volume	0.0182 cbm
Weight of case	3.64 kg
EAN of inner box	4038917925174
EAN case	4038917925181

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

