



Product description

| | |
|---------------------|---|
| Order number | 65.716.302 |
| Product description | Screw cap, yellow, suitable for screw cap micro tubes, 500 piece(s)/bag |

Product characteristics

| | |
|--------------|-----------------------|
| Cap | flat cap |
| Closure type | screw cap |
| Suitable for | screw cap micro tubes |

Material & colours

| | |
|-------------------|--|
| Product material | Polypropylene (PP) |
| Colour of product | yellow |
| O-ring material | Ethylene propylene diene rubber (EPDM) |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap, yellow, suitable for screw cap micro tubes

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 10000 |
| Piece(s) / pallet | 720000 |
| EAN of inner box | 4038917925549 |
| EAN case | 4038917925556 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

