



**Product Specification**  
**Screw cap, violet, suitable for screw cap micro tubes**

**Page 1**



**Product description**

|                     |   |
|---------------------|---|
| Order number        | 65.716.308  |
| Product description | Screw cap, violet, suitable for screw cap micro tubes, 500 piece(s)/bag |

**Product characteristics**

|              |                       |
|--------------|-----------------------|
| Cap          | flat cap              |
| Closure type | screw cap             |
| Suitable for | screw cap micro tubes |

**Material & colours**

|                   |  |
|-------------------|--|
| Product material  | Polypropylene (PP)                     |
| Colour of product | violet                                 |
| O-ring material   | Ethylene propylene diene rubber (EPDM) |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap, violet, suitable for screw cap micro tubes**

**Page 2**

**Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Batched          | yes                           |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 10000         |
| Piece(s) / pallet             | 720000        |
| EAN of inner box              | 4038917925389 |
| EAN case                      | 4038917925396 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**