



Product Specification
Screw cap, red, sterile, suitable for screw cap micro tubes

Page 1



Product description

| | |
|---------------------|--|
| Order number | 65.716.323 |
| Product description | Screw cap, red, sterile, suitable for screw cap micro tubes, 500 piece(s)/double bag |

Product characteristics

| | |
|--------------|-----------------------|
| Cap | flat cap |
| Closure type | screw cap |
| Suitable for | screw cap micro tubes |

Material & colours

| | |
|-------------------|--|
| Product material | Polypropylene (PP) |
| Colour of product | red |
| O-ring material | Ethylene propylene diene rubber (EPDM) |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-06

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Screw cap, red, sterile, suitable for screw cap
micro tubes**

Page 2

Purity & certification

| | |
|------------------|--------------------------------|
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Type of smallest subpackaging | double bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 5000 |
| Piece(s) / pallet | 540000 |
| Depth of case | 351 mm |
| Width of case | 280 mm |
| Height of case | 138 mm |
| Case volume | 0.0136 cbm |
| Weight of case | 2.05 kg |
| EAN of inner box | 4038917923415 |
| EAN case | 4038917923422 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-06

This document was prepared by EDP support and is valid without signature.



SARSTEDT