



**Product Specification**  
**Screw cap, black, sterile, suitable for screw cap micro tubes**

**Page 1**



**Product description**

Order number	65.716.728
Product description	Screw cap, black, sterile, suitable for screw cap micro tubes, 500 piece(s)/double bag

**Product characteristics**

Cap	flat cap
Closure type	screw cap
Suitable for	screw cap micro tubes

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	black
O-ring material	Ethylene propylene diene rubber (EPDM)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap, black, sterile, suitable for screw cap  
micro tubes**

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	double bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	390000
Depth of case	418 mm
Width of case	278 mm
Height of case	130 mm
Case volume	0.0151 cbm
Weight of product	0.0003 kg
Weight of case	2.01 kg
EAN of inner box	4038917119016
EAN case	4038917053624

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

