



Product description

| | |
|---------------------|--|
| Order number | 72.690.001 |
| Product description | Reaction tube, working volume: 1.5 ml, material: PP, transparent, cap: transparent, retaining loop cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------|--------------------|
| Graduation | yes |
| Cap | retaining loop cap |
| Volume of work | 1,5 ml |
| Centrifugation max (RCF) | 20000 x g |

Size

| | |
|-------------------|---------|
| Diameter | 10.8 mm |
| Length of product | 39 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | Polypropylene (PP) |
| Colour of cap | transparent |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 5000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 5000 |
| Piece(s) / pallet | 150000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 340 mm |
| Case volume | 0.0539 cbm |
| Weight of product | 0.001 kg |
| Weight of case | 5.56 kg |
| EAN of inner box | 4038917161657 |
| EAN case | 4038917554138 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

