



Product description

Order number	72.690.307
Product description	Reaction tube, working volume: 1.5 ml, material: PP, purple, cap: purple, cap attached, with injected graduation and writing space, 500 piece(s)/bag

Product characteristics

Graduation	yes
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g
Autoclavable (max. °C)	121 °C

Size

Sample volume	1.5 ml
Length of product	39 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	purple
Colour of cap	purple

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	160000
Case volume	0.0558 cbm
Weight of product	0.001 kg
Weight of case	5.6 kg
EAN of inner box	4038917369329
EAN case	4038917369312

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

