



### **Product description**

Order number	72.691
Product description	Reaction tube, working volume: 2 ml, material: PP, transparent, cap: transparent, retaining loop cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag

### **Product characteristics**

Graduation	yes
Cap	retaining loop cap
Closure type	hinged lid
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

### **Size**

Diameter	10.8 mm
Length of product	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	transparent

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	120000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0011 kg
Weight of case	6.05 kg
EAN of inner box	4038917312790
EAN case	4038917312783

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

