



Product description

Order number	72.698.200
Product description	Reaction tube, working volume: 0.5 ml, material: PP, transparent, without cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	7.8 mm
Length of product	30 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	330000
Depth of case	438 mm
Width of case	268 mm
Height of case	155 mm
Case volume	0.0182 cbm
Weight of product	0.0004 kg
Weight of case	2.3 kg
EAN of inner box	4038917118118
EAN case	4038917015875

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

