



### **Product description**

Order number	72.704.001
Product description	SafeSeal reaction tube, working volume: 0.5 ml, material: PP, brown, cap: brown, SafeSeal cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag

### **Product characteristics**

Graduation	yes
Cap	SafeSeal cap
Volume of work	0,5 ml
Centrifugation max (RCF)	30000 x g

### **Size**

Diameter	7.8 mm
Length of product	30 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	brown
Closure material	Polypropylene (PP)
Colour of cap	brown

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	156000
Depth of box	189 mm
Width of box	151 mm
Height of box	97 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0005 kg
Weight of case	1.61 kg
EAN of inner box	4038917227520
EAN case	4038917227513

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

