



Product description

| | |
|---------------------|---|
| Order number | 72.704.006 |
| Product description | SafeSeal reaction tube, working volume: 0.5 ml, material: PP, blue, cap: blue, SafeSeal cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------|--------------|
| Graduation | yes |
| Cap | SafeSeal cap |
| Volume of work | 0,5 ml |
| Centrifugation max (RCF) | 30000 x g |

Size

| | |
|-------------------|--------|
| Diameter | 7.8 mm |
| Length of product | 30 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | blue |
| Closure material | Polypropylene (PP) |
| Colour of cap | blue |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 156000 |
| Depth of box | 189 mm |
| Width of box | 151 mm |
| Height of box | 97 mm |
| Depth of case | 388 mm |
| Width of case | 358 mm |
| Height of case | 135 mm |
| Case volume | 0.0188 cbm |
| Weight of product | 0.0005 kg |
| Weight of case | 1.61 kg |
| EAN of inner box | 4038917227643 |
| EAN case | 4038917227636 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

