



Product description

Order number	72.609.300
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, without cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Knurls	yes
Base shape	skirted conical base
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Satisfies the requirement	IATA
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	180000
Case volume	0.0492 cbm
Weight of product	0.0012 kg
Weight of case	6.83 kg
EAN of inner box	4038917434287
EAN case	4038917434270

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

