



### Product description

Order number	72.694.323
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, cap: red, cap assembled, with printed writing space, with graduations, sterile, 100 piece(s)/bag

### Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Knurls	yes
Base shape	skirted conical base
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Size**

Diameter	10.8 mm
Length of product	44 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	red
O-ring material	Ethylene propylene diene rubber (EPDM)

### **Purity & certification**

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### **Packaging**

Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	108000
Depth of case	351 mm
Width of case	280 mm
Height of case	138 mm
Case volume	0.0136 cbm
Weight of case	1.82 kg
EAN of inner box	4038917926577
EAN case	4038917926584

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

