



Product description

Order number	72.730.007
Product description	Screw cap micro tube, working volume: 0.5 ml, skirted conical base, with knurling, transparent, cap: natural, cap enclosed, with printed writing space, 250 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	no
Cap	standard screw cap
Knurls	yes
Base shape	skirted conical base
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Diameter	10.8 mm
Length of product	44 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

Purity & certification

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	5000
Piece(s) / pallet	150000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0015 kg
Weight of case	8.18 kg
EAN of inner box	4038917119290
EAN case	4038917016568

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

