




## Product description

|                     |                                                                                                                                                                                                                                                                                                           |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 93.873.681                                                                                                                                                                                                                                                                                                |
| Product description | Cryobox, with numerical coding at each aperture, for low-temperature storage, material: PC, blue, slip-on lid with ventilation function, cap: transparent, (LxWxH): 132 x 132 x 53 mm, format: 9 x 9, for 81 collection tubes, for CryoPure tubes 1.2?2.0 ml internal and external thread, 5 piece(s)/bag |

## Product characteristics

|                  |                                                        |
|------------------|--------------------------------------------------------|
| Graduation       | no                                                     |
| Cap              | slip-on lid with ventilation function                  |
| Closure type     | slip-on lid                                            |
| Design           | with numerical coding at each aperture                 |
| Suitable for     | CryoPure tubes 1.2?2.0 ml internal and external thread |
| Storage capacity | 81                                                     |
| Format           | 9 x 9                                                  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

|                   |        |
|-------------------|--------|
| Width of product  | 132 mm |
| Height of product | 53 mm  |
| Length of product | 132 mm |

### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polycarbonate (PC) |
| Colour of product | blue               |
| Colour of cap     | transparent        |

### **Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 5             |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 5             |
| Piece (s) / outer case        | 20            |
| Piece(s) / pallet             | 720           |
| Depth of case                 | 135 mm        |
| Width of case                 | 270 mm        |
| Height of case                | 540 mm        |
| Case volume                   | 0.0197 cbm    |
| Weight of product             | 0.1694 kg     |
| Weight of case                | 3.91 kg       |
| EAN of inner box              | 4038917122887 |
| EAN case                      | 4038917415064 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

