




**Product Specification**  
**Storage box, hinged lid, PP, format: 10 x 10, for 100 collection tubes**

**Page 1**



**Product description**

Order number	93.877.510
Product description	Storage box, hinged lid, material: PP, green, format: 10 x 10, for 100 collection tubes, suitable for vessels with maximum dimensions of 45 x 12 mm, 5 piece(s)/bag

**Product characteristics**

Closure type	hinged lid
Suitable for	vessels with maximum dimensions of 45 x 12 mm
Storage capacity	100
Format	10 x 10
Autoclavable (max. °C)	121 °C

**Size**

Width of product	140 mm
Height of product	55 mm
Length of product	150 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	green
Closure material	Polypropylene (PP)
Colour of cap	green

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	20
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	20
Piece(s) / pallet	720
Depth of case	305 mm
Width of case	279 mm
Height of case	305 mm
Case volume	0.026 cbm
Weight of product	0.1835 kg
Weight of case	4.1 kg
EAN of inner box	4038917161237
EAN case	4038917161220

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

