



Product description

| | |
|---------------------|--|
| Order number | 93.877.999 |
| Product description | Storage box, hinged lid, material: PP, colour mix, format: 10 x 10, for 100 collection tubes, suitable for vessels with maximum dimensions of 45 x 12 mm |

Product characteristics

| | |
|------------------------|---|
| Closure type | hinged lid |
| Suitable for | vessels with maximum dimensions of 45 x 12 mm |
| Storage capacity | 100 |
| Format | 10 x 10 |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|-------------------|--------|
| Width of product | 140 mm |
| Height of product | 55 mm |
| Length of product | 150 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | colour mix |
| Closure material | Polypropylene (PP) |
| Colour of cap | colour mix |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|------------------------|---------------|
| Minimum order qty. | 5 |
| Piece(s) / inner box | 5 |
| Piece (s) / outer case | 20 |
| Piece(s) / pallet | 720 |
| Case volume | 0.026 cbm |
| Weight of product | 0.204 kg |
| Weight of case | 4.09 kg |
| EAN of inner box | 4038917240376 |
| EAN case | 4038917240369 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

