



Product description

| | |
|---------------------|--|
| Order number | 93.879 |
| Product description | Rack, material: PP, multi-coloured, format: 12 x 8, (LxWxH): 108 x 108 x 257 mm, suitable for micro tubes and tubes with Ø 4-30 mm, 1 piece(s)/bag |

Product characteristics

| | |
|------------------------|--------------------------------------|
| Suitable for | micro tubes and tubes with Ø 4-30 mm |
| Format | 12 x 8 |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|-------------------|--------|
| Width of product | 108 mm |
| Height of product | 257 mm |
| Length of product | 108 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Rack, PP, format: 12 x 8, suitable for micro tubes
and tubes with Ø 4-30 mm

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | multi-coloured |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 4 |
| Piece(s) / pallet | 384 |
| Depth of case | 216 mm |
| Width of case | 254 mm |
| Height of case | 241 mm |
| Case volume | 0.0132 cbm |
| Weight of product | 0.8348 kg |
| Weight of case | 3.6 kg |
| EAN of inner box | 4038917355162 |
| EAN case | 4038917355155 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

