



### **Product description**

Order number	95.064.259
Product description	Box, Styrofoam, 20-place, for micro tubes, 8 pcs./pack

### **Product characteristics**

Storage capacity	20
Format	5 x 4

### **Material & colours**

Product material	Styrofoam
Colour of product	white

### **Purity & certification**

Product category	No marking necessary
Batched	Lot number

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	16
Piece(s) / inner box	8
Piece (s) / outer case	56
Piece(s) / pallet	1792
Depth of case	401 mm
Width of case	451 mm
Height of case	280 mm
Case volume	0.0504 cbm
Weight of case	1.36 kg
EAN of inner box	4038917386883
EAN case	4038917386876

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

