



Product description

| | |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number | 95.64.922 |
| Product description | Storage box, slip-on lid, material: cardboard, white, format: 7 x 7, for 49 collection tubes, suitable for tubes up to max. 22 mm Ø, tube height between 88.5-129 mm, 1 piece(s)/bag |

Product characteristics

| | |
|------------------|-----------------------------------------------------------|
| Closure type | slip-on lid |
| Suitable for | tubes up to max. 22 mm Ø, tube height between 88.5-129 mm |
| Storage capacity | 49 |
| Format | 7 x 7 |

Size

| | |
|-------------------|--------|
| Width of product | 154 mm |
| Height of product | 130 mm |
| Length of product | 152 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|-----------|
| Product material | cardboard |
| Colour of product | white |
| Closure material | cardboard |
| Colour of cap | white |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 27 |
| Piece(s) / pallet | 324 |
| Depth of case | 521 mm |
| Width of case | 521 mm |
| Height of case | 610 mm |
| Case volume | 0.1653 cbm |
| Weight of product | 0.22 kg |
| Weight of case | 7.3 kg |
| EAN of inner box | 4038917095709 |
| EAN case | 4038917067447 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

